


# HILBRO

The Surgical Instruments Manufacturing Company




# Omni Tract

RETRACTOR SYSTEM (ONE PIECE)

Omni-Tract Surgical has over 30 years of experience in table mounted retractor design.

The concept of placing the retractor frame to the O.R. table rail allows you to place retractor blades at any angle with exact pressure for unilateral, bilateral, single or multi-level exposures. Table mounted retraction eliminates hand held retractors and self retaining retractor migration. Wishbone retractor systems fit all patient sizes, and can be used in multiple surgical specialties.



**Description:**

**What's Included:**

- Pan/Cover Ser; Qty: 1
- Balfour Center Blade, Horgan, 15 3/4" (400mm), Malleable, 2" Wide x 7 3/4" Deep; Qty: 1
- Abdominal Retractor 18-1/2"; Qty: 1
- Balfour Fourth Blade Attachment, For Right Arm, Arm Depth 241mm; Qty: 1
- Balfour Center Blade; Qty: 1
- Balfour Wire Lateral Blades; Qty: 4
- Gibson Balfour Retractor; Qty: 1
- Flex Clamp; Qty: 1
- Lapo-Clamp; (REF 3720) Qty: 1
- Clamp Sterile Field; (REF 3700) Qty: 1
- Retractor Frame Component, Angled Arm, 55 Degree, 7"x16"; (REF 3732) Qty: 1

**HILBRO**

The Surgical Instruments Manufacturing Company

Small Industries Estate, Sialkot-51310, Pakistan

Phone: +92 (52) 3550922, 3255459

Fax: +92 (52) 3553439, 6526579

e-mail: sales@hilbro-instruments.com

web: hilbro.net